

Spinsoft

Care at its very best!



Spinal Needle

Quincke & Epidot Pencil Point


Spinsoft

Zhejiang FERT mainly manufactures disposable infusion set with precise filter and anesthesia puncture kits and other similar products.

Zhejiang FERT, was honored as TongXiang High-Tech enterprise in 2002. Jiaxing High-Tech enterprise in 2003 and Zhejiang High-Tech enterprise in 2005, FERT also took the lead in passing the inspection by S. F. D. A. (State Food and Drug Administration in China) for its disposable Anesthesia Manufacture Implementation in November 2003, and was certified ISO 9001: 2000, ISO 13485: 2003 and obtained its first CE marking in June 2004.

QUINCKE

- Smooth but moderate insertion resistance of the needle makes it possible to identify easily the tissue site through where the needle is inserted.
- The proven trihedral tip combines smooth entry with just the right amount of puncture resistance to assist safe dural puncture.
- The base of the needle is transparent, making it easy to check for reflux of spinal fluid.



EPIDOT PENCIL POINT

- The conical design of the needle point makes it possible to minimize lesions of the dura fibers and to reduce the incidence of postoperative headaches.



EPIDURAL ANESTHESIA NEEDLES “FIXED WING”

Features:

- The ultrathin wall of cannula and the curved design of the needle point allows protection of the catheter from damage during procedures and smooth insertion into the epidural space.



QUINCKE AND EPIDOT PENCIL POINT

GAUGE CALIBRE	LENGTH (mm)	COLOUR	GUIDE NEEDLE (METAL OR PLASTIC HUB)
25 G	70-90mm	ORANGE	21 G X 38mm
27 G	70-90mm	GRAY	22 G X 38mm

Marketed by:



+92-21-345 411 17-8